

UV-Vis Spectrophotometer

UV-1900i Plus



SHIMADZU

Log Out Mode Menu Ready 12/24 16:49

 Photometric	 Spectrum	 Quantitation
 Kinetics	 Time Course	 Bio Method

Method File Data File Method DC Control Control Software

UV-1900i Plus UV-VIS SPECTROPHOTOMETER

Navigate Your Way

Refined User-Friendliness

Easy-to-use touch-screen display

Faster response speed with adoption of a new CPU

Various Functions for Comfortable Daily Measurement

User support with assist and shutdown/wakeup functions

High Performance to Meet Diverse Needs

The instrument is equipped with an ultra high-speed scan, which can acquire a spectrum in a few seconds, with the lowest level of stray light and noise in its class

Validation functions enable checks in accordance with Pharmacopoeia (JP/USP/EP) to be performed easily*

* Operations are performed via the touch panel of the UV-1900i Plus. If performing validation in compliance with USP/EP using PC software, optional validation software is required.

UV-i Selection



UV-1900i Plus



UV-2600i Plus
UV-2700i Plus



UV-3600i Plus



SolidSpec™-3700i

Refined User-Friendliness




The instrument is equipped with a stylus pen. Operation can be performed with this pen or a finger.

Easy-to-Use Interface

Grasp the Current Status and Operating Procedures at a Glance

The UV-1900i Plus on-screen user interface includes large, easy-to-see icons deployed on a black background, so the instrument settings are evident at a glance. In addition, the large, easy-to-see icons improve intuitive understanding, which enables users to quickly become familiar with the operations. Furthermore, the user interface is designed to minimize transitions between windows, so users do not get confused during the operations. The system also features a faster response speed by adopting a new CPU.



Display languages are available in eight languages (Japanese, English, Chinese, Spanish, Portuguese, German, French, Russian).

Multiple Measurement Modes

UV-1900i Plus has six basic measurement modes, which offer the optimal solution for a variety of measurements.


Photometric

Measures the photometric value at a single wavelength or multiple (up to eight) wavelengths.








Spectrum

Measures a sample spectrum using wavelength scanning.



Quantitation

Generates a calibration curve from the measurement of standards, and then calculates the concentrations of unknowns.



Kinetics

Measures absorbance changes as a function of time, and obtains the enzymatic activity values. The kinetics measurement method or the rate measurement method can be selected.

Time Course

Measures changes over time in photometric values at a specified wavelength.

Biomethod

Quantifies DNA or protein concentrations.

* The background color of the operation screen can also be set to white.

Various Functions for Comfortable Daily Measurement

Assist Function

The instrument assists the user to ensure that measurements are performed using the correct procedure.*1



Checks whether or not baseline correction, auto zero correction or cell blank correction has been performed and informs the user if it has not been.

Informs the user if the last correction performed is not appropriate for the planned measurement.


Informs the user if the UV-1900i Plus is not ready to begin measurement and when starting measurements and 100 %T (0 Abs) corrections.

*1 This function can be enabled/disabled.

Shutdown Function


The shutdown function puts UV-1900i Plus to sleep after it's been used, or after a certain period of time. Putting the instrument into sleep mode limits power consumption and helps preserve the lamp.

Configured to shut down after 2 hours



Wakeup Function

Automatically wakes up a sleeping UV-1900i Plus at a specified time. This function eliminates the need for analysts to wait for it to warm up, allowing them to be more productive.



- Automated support functions utilizing digital technology, such as M2M, IoT, and Artificial Intelligence (AI), that enable higher productivity and maximum reliability.
- Allows a system to monitor and diagnose itself, handle any issues during data acquisition without user input, and automatically behave as if it were operated by an expert.
- Supports the acquisition of high quality, reproducible data regardless of an operator's skill level for both routine and demanding applications.

Startup Validation Function

Instrument performance checks can be performed automatically when the instrument is turned ON.*2

This enables even more reliable instrument operation.

It can also be implemented in combination with the Wakeup function.


*2 This function can be enabled/disabled.

The performance check items that can be used with this function are as follows: wavelength accuracy (D2), wavelength repeat accuracy (D2), degradation, noise level, baseline flatness, baseline stability (drift), [EP] Control of Wavelength Accuracy (D2), and [USP] Control of Wavelengths (D2).



Network Connectivity Function

Data can now be transferred to a PC via a network. With wireless printing, multiple UV units can print from a single printer.*3 (A router and other network equipment must be installed to use a network.)



*3 A router with wireless function is required.

*4 Optional expanded memory is required. The instrument is not compatible with control via a network.

*5 A PictBridge compatible printer is required.


Bar-Code Reader and Keyboard Entry Function

Sample names and numerical values can be entered by a bar-code reader or from the keyboard.

This saves time when entering sample names for a multi-sample analysis, and prevents sample misidentification and other human errors.*6



*6 Use a bar-code reader and keyboard with a USB connection.




High Performance to Meet Diverse Needs



Low Stray Light

Linearity will be lost in the high absorbance region due to stray light. UV-1900 series is equipped with Low Stray Light Diffraction Gratings, achieve "best-in-class" low stray light. Stray light is at 0.5 % max. (198 nm), making accurate measurements possible at 1 Abs or more, even in the ultraviolet region. In addition, high-concentration samples can be quantified accurately.


The figure on the right is a calibration curve for acetic acid, created with absorbance at 200 nm. The correlation coefficient is 0.9998 and correct measured values are obtained even with high-concentration samples.




Ultra-Fast Scan

Spectra can be acquired as fast as 29,000 nm/min. Ultra-fast scanning is effective for tracking chemical reactions in a short time. In addition to the absorbance change at specified wavelengths, spectra can be acquired in a short time with the UV-1900i Plus. Therefore, more detailed behavior can be investigated by observing spectra with the UV-1900i Plus.

The figures below show the analysis of the particle agglomeration process when salts are added to silver nanoparticles. Measurements of the 300 to 700 nm region were performed in ultra-fast scan mode. In addition to the decrease of absorbance at 400 nm and the increase of absorbance at 480 nm, the temporal changes of spectra can be observed.



Absorption Spectra of Silver Nanoparticles




Temporal Changes of Absorbance at 400 and 480 nm

High Reproducibility and Repeatability Accuracy

The photometric repeatability accuracy is 0.0002 Abs max. (0.5 Abs and 1.0 Abs). As a result, variance in the measurement results is suppressed, enabling more accurate quantitation and the detection of low-concentration samples.

The figure on the right is a calibration curve for caffeine, created with absorbance at 273 nm. The calibration curve has an Abs = 0.0528 Conc. The lower limit of quantitation determined from the standard deviation is 0.0051 mg/L.

No.	Absorbance of Blank Solution (273 nm)
1	-0.00001
2	0.00001
3	-0.00002
4	0.00002
5	0.00001
6	-0.00003
7	0.00001
8	-0.00004
9	0.00001
10	0.00005
Standard Deviation σ	0.000025



Note: One method of determining the lower limit of quantitation is to use ten times the standard deviation. This is an actual measured value and is not guaranteed.

Enables Compliance with ER/ES Regulations and Stronger Data Management

LabSolutions™ UV-Vis Software

Enables higher productivity and provides for a more convenient analytical environment.



Startup

Different measurement applications can be started from the application launcher. Extensive UV validation software programs supporting instrument performance checks as well as optional software programs for special analyses are also available, so various measurements can be intuitively performed.




Four Measurement Mode Windows

Setting Parameters

Instrument Control Panel

Instrument parameter settings can be specified via panels that are separate from the measurement window. The control panels include various functionality that is laid out for superior visibility. Each measurement window connects seamlessly to the corresponding parameter settings window.



Presets (Spectrum and Time Course)

When multiple samples are measured, the sample information can be configured in advance.

Sample Information	Type	Sample Name	Sample ID	File Path	File Name	Parameter	Comments
1	Sample	Sample001	Sample001	Path\Sample001	Sample001	UV-Vis	
2	Sample	Sample002	Sample002	Path\Sample002	Sample002	UV-Vis	
3	Sample	Sample003	Sample003	Path\Sample003	Sample003	UV-Vis	
4	Sample	Sample004	Sample004	Path\Sample004	Sample004	UV-Vis	
5	Sample	Sample005	Sample005	Path\Sample005	Sample005	UV-Vis	
6	Sample	Sample006	Sample006	Path\Sample006	Sample006	UV-Vis	
7	Sample	Sample007	Sample007	Path\Sample007	Sample007	UV-Vis	
8	Sample	Sample008	Sample008	Path\Sample008	Sample008	UV-Vis	
9	Sample	Sample009	Sample009	Path\Sample009	Sample009	UV-Vis	
10	Sample	Sample010	Sample010	Path\Sample010	Sample010	UV-Vis	
11	Sample	Sample011	Sample011	Path\Sample011	Sample011	UV-Vis	
12	Sample	Sample012	Sample012	Path\Sample012	Sample012	UV-Vis	
13	Sample	Sample013	Sample013	Path\Sample013	Sample013	UV-Vis	
14	Sample	Sample014	Sample014	Path\Sample014	Sample014	UV-Vis	
15	Sample	Sample015	Sample015	Path\Sample015	Sample015	UV-Vis	
16	Sample	Sample016	Sample016	Path\Sample016	Sample016	UV-Vis	
17	Sample	Sample017	Sample017	Path\Sample017	Sample017	UV-Vis	
18	Sample	Sample018	Sample018	Path\Sample018	Sample018	UV-Vis	
19	Sample	Sample019	Sample019	Path\Sample019	Sample019	UV-Vis	
20	Sample	Sample020	Sample020	Path\Sample020	Sample020	UV-Vis	

Instrument Check

Assist Function* ANALYTICAL INTELLIGENCE

The user can tell at a glance if the instrument has finished warming up. Also, the user is notified when attempting to implement measurements without making necessary corrections, thereby supporting appropriate analysis work.

* This function can be enabled/disabled.



Reliable LabSolutions Software

In addition to LabSolutions UV-Vis, which provides basic functionality, Shimadzu offers LabSolutions DB UV-Vis and LabSolutions CS UV-Vis to meet the requirements of ER/ES regulations.


LabSolutions DB UV-Vis

LabSolutions DB UV-Vis allows for secure data management by integrating a data management function with LabSolutions UV-Vis. Compliant with ER/ES regulations, the software is optimally configured for customers using a PC. It is recommended for facilities that do not require network connections and want to be ER/ES compliant.



LabSolutions CS UV-Vis


LabSolutions CS, which is freely accessible to the analysis network, can be connected to LabSolutions UV-Vis, eliminating the need for connecting a PC to the instrument. Since all the data are managed on a server, LabSolutions CS UV-Vis can be read from any personal computer on a network. With terminal service, LabSolutions UV-Vis can be controlled from a client PC without installing LabSolutions UV-Vis on it. It is recommended for facilities that have a large number of users, manage data in a database, and want to be ER/ES compliant.




^{*1} The acquisition controller PC controls analytical instruments.
^{*2} A terminal server is a server for using terminal services. Users can view data reports and perform electronic signature operations through terminal services. It is ideal for remote connections because of the low network load. Only LC, GC, LCMS, and GCMS support analysis and post-run operations through terminal services.
^{*3} If a terminal service is used, LabSolutions software does not need to be installed on client PCs or tablets.
^{*4} Servers can be built on various clouds (IaaS). AWS (Amazon Web Services), Microsoft® Azure®, GCP™ (Google Cloud Platform™)

Validation Functions Compliant with JP, USP, and EP (Touch-screen Display)

This instrument can not only run checks for nine JIS items, but also those stipulated in the Japanese Pharmacopoeia (JP), United States Pharmacopoeia (USP), and European Pharmacopoeia (EP). The UV-1900 series also complies with the hardware specifications required by each Pharmacopoeia. In addition, the check conditions can be saved. As a result, checks can be performed easily just by calling them up as needed. Check results can also be saved.



Main Screen




Parameters Setting Screen



Check Results Screen

Improved Security Functions

An external control security function has been added to provide more support for compliance with regulations. Three user authority levels, "Administrator", "Developer", and "Operator", can be set for instrument users.





Resolution of 1 nm, the Highest in its Class

In addition to achieving a resolution of 1 nm, the highest in its class, by using a monochromator with a Czerny-Turner mounting, the UV-1900i Plus also features a compact, bright optical system. The instrument is more than capable of meeting the wavelength resolution required in the USP and EP.

Compatible with Validation from PC Software

UV validation software (limited function edition) is included as standard. Using this software, inspections can be implemented as stipulated for 9 items in JIS and the JP. If optional UV validation software (USP/EP compatible edition) is added, examinations can be implemented as stipulated in the USP and EP.



- Inspection results can not only be printed but also saved to a file, so the results can be called up later for confirmation.
- The inspection parameters can also be saved to separate files for periodic and routine inspections, and then called up for use.


- The user can select confirmation of instrument performance indicators as per JIS K0115 General rules for molecular absorptiometric analysis, as well as the general test methods in the JP. (Order inspection jigs and reagents separately.)

Note: Optional validation software is required to support the USP and EP.

Applications



Foods

This is an example of the analysis of food dyes. By using ultra-fast scan mode, the time needed for routine spectral checks can be shortened. The 350 nm to 900 nm region can be measured at 1 nm intervals in approx. 4 seconds.





Pharmaceuticals and Life Sciences

This is an example of the analysis of λ DNA. Trace quantities (on the order of a few μ L) can be measured by combining the instrument with Nano Stick and TrayCell[®]. Using TrayCell, a calibration curve for 4 μ L of λ DNA was obtained in the range between 25 ng/ μ L and 500 ng/ μ L.



Chemistry

Materials can be confirmed quantitatively by using LabSolutions UV-Vis and color measurement software. LabSolutions UV-Vis color software shows the chromaticity diagram of the XY color system or Hazen color.



Accessories

Expanded Memory (P/N 207-23119-41)

The UV-1900i Plus main unit can store up to 999 sets of data. The saved data can be read out from a PC. (A network connection is required for data readout.)

Film Holder (P/N 204-58909)

Used in transmittance measurement of thin samples such as films and filters. It holds the samples for analysis.



Long-Path Rectangular Cell Holder (P/N 204-23118-01)

Holds two rectangular cells with an optical path length of 10, 20, 30, 50, 70, or 100 mm.



Test Tube Holder (P/N 207-23510-41)

Holds test tube instead of sample compartment.

- Outside diameter: \varnothing 15 to 22 mm
- Height: 75 to 115 mm

Note: Test tube is not included.



CPS-100 Cell Positioner, Thermoelectrically Temperature Controlled (P/N 206-29500-42, -43, -58)

This attachment permits measurement of up to six sample cells under constant-temperature conditions. Combination of this attachment and the Kinetics mode provides measurement of temperature-sensitive enzyme kinetics of one to six samples.

- Number of cells: 6 on the sample side (temperature-controlled)
1 on the reference side (temperature not controlled)
- Temperature control range: 16 to 60 °C
- Temperature display accuracy (difference from the true value): ± 0.5 °C
- Temperature control precision (variation of temperature): ± 0.1 °C
- Ambient temperature: 15 to 35 °C

Note: Square cells (P/N 200-34442) are not included, please purchase separately. A USB adapter CPS (P/N 206-25234-91) is required



Multi-Cell Sample Compartment (Type-p) (P/N 206-69160-58)

Holds up to six 10-mm square cells on the sample side. No temperature control capability.




- Number of cells: 6 on the sample side
1 on the reference side

Note: Square cells are not included, please purchase separately.

Sipper Unit

Model	P/N	Standard Sample Volume
160L (Standard Sipper)	206-23790-51	2.0 mL
160T (Triple-Pass Sipper)	206-23790-52	1.5 mL
160C (Constant-Temperature Sipper)	206-23790-53	2.5 mL
160U (Supermicro Sipper)	206-23790-54	0.5 mL

Four types of sipper units with different flow cells are available. The stepping motordriven peristaltic pump ensures reliable and smooth aspiration of sample solution.



(Direct driving is possible from the UV-1900i Plus so no interface is required.)

Note: The use of a Solenoid Valve (Fluoropolymer) (P/N 204-06599-01) and the SWA-2 Sample Waste Unit (P/N 206-23820-58) are recommended when strong acids, strong alkalis, or organic solvents are to be measured.

TCC-100 Thermoelectrically Temperature Controlled Cell Holder (P/N 206-29510-42, -43, -58)


Uses Peltier effect for controlling the sample and reference temperature, so no thermostated bath or cooling water is required.

- Number of cells: One each on the sample and reference sides (temperature-controlled)
- Temperature control range: 7 to 60 °C
- Temperature display accuracy (difference from the true value): ± 0.5 °C
- Temperature control precision (variation of temperature): ± 0.1 °C

Note: Square cells (P/N 200-34442) are not included, please purchase separately.



Various other accessories



SolidSpec, the Analytical Intelligence Logo, LabSolutions and PPSQ are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries. Excel, Microsoft and Azure are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Amazon Web Services and AWS are trademarks of Amazon.com, Inc. or its subsidiaries. Google Cloud Platform and GCP are trademarks of Google LLC. TrayCell is a registered trademark of Hellma GmbH & Co. KG. in Germany.



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.